









# DIAMOND®-Tips








A large choice of packaging to fit your activity and budget

Diamond certified tips assure optimal performance of your Pipetman pipette.



Tip model	Pipette models U = Pipetman Ultra, C = Pipetman Concept	Volume range	Tipack Diamond Tips:		Tower Pack: refill system for Diamond Tips		SterilPack: 400 tips individually wrapped	EcoPack: box of 10,000 tips (10 bags of 1,000)	EasyPack: box of 1,000 tips (5 bags of 200)	Empty Racks (set of 2)				
			10 Tipacks of Autoclavable	96 racked tips $\gamma$ -Ray sterilized	10 racks of 96 tips Autoclavable	10 racks of 96 tips $\gamma$ -Ray sterilized								
 <b>D10</b>	P2, P10, U2, U10, C10	0.2 - 10 $\mu$ L	F161632	F161905	F167101	F167201	F161310	F161630	F161631	F161791				
 <b>DL10</b>	P2, P10, U2, U10, U8x20, U12x20, C10, C8x10, C12x10	0.2 - 20 $\mu$ L	F161452	F161455	F167102	F167202	F161320	F161450	F161451	F161393				
 <b>D200</b>	P20 to P200, P8x200, F2 to F200, U20 to U200, U8x20, U12x20, U8x300, U12x300, C100, C300, C8x100, C12x100, C8x300, C12x300	2 - 200 $\mu$ L	F161932	F161937	F167103	F167203	F161330	F161930	F161931	F161391				
 <b>D300</b>	P200, U200, U8x300, U12x300, C300, C8x300, C12x300	20 - 300 $\mu$ L	F161732	F161735	F167105	F167205	-	F161730	F161731	F161737				
 <b>D1000</b>	P1000, F250 to F1000, U1000	200 - 1000 $\mu$ L	F161672	F161975	F167104 <small>(7 racks of 96 tips)</small>	F167204 <small>(7 racks of 96 tips)</small>	F161340	F161670	F161671	F161793				
 <b>D1200</b>	P1000, F250 to F1000, U1000, C1200, C8x1200, C12x1200	100 - 1200 $\mu$ L	F161102	F161103	-	-	-	F161100	F161101	F106104				
 <b>D5000</b>	P5000, U5000, C5000	1 - 5 mL	F161370 <small>(12 Tipacks of 50 tips)</small>	-	<table border="1"> <thead> <tr> <th>Item</th> <th>Reference Number</th> </tr> </thead> <tbody> <tr> <td>Universal reload box</td> <td>F167100</td> </tr> </tbody> </table>		Item	Reference Number	Universal reload box	F167100	-	-	F161571 <small>(3 bags of 334)</small>	F161360 <small>(1 rack)</small>
Item	Reference Number													
Universal reload box	F167100													
 <b>D10ml</b>	P10ml, U10ml, C10ml	1 - 10 mL	F161230 <small>(8 Tipacks of 30 tips)</small>	-	-	-	-	-	F161210 <small>(200 tips)</small>	F161220 <small>(1 rack)</small>				


## Work safety with the Diamond Filter Tips

 <b>DF10</b>	P2, P10, U2, U10	0.2 - 10 $\mu$ L	F161710	F161633
 <b>DFL10</b>	P2, P10, U2, U10, U8x20, U12x20, C10, C8x10, C12x10	0.2 - 10 $\mu$ L	F161454	F161453
 <b>DF30</b>	P20, F2 to F25, U20, U8x20, U12x20	2 - 30 $\mu$ L	F161935	F161933
 <b>DF100</b>	P100, F50 to F100, U100, C100, C8x100, C12x100	20 - 100 $\mu$ L	F161939	F161938
 <b>DF200</b>	P200, P8x200, F200, U200, U8x300, U12x300, C300, C8x300, C12x300	20 - 200 $\mu$ L	F161936	F161934
 <b>DF300</b>	P200, U200, U8x300, U12x300, C300, C8x300, C12x300	20 - 300 $\mu$ L	F161734	F161733
 <b>DF1000</b>	P1000, F250 to F1000, U1000	200 - 1000 $\mu$ L	F161770	F161673

The only way to ensure precision and reproducibility across the full volume range is to use Diamond precision tips with Pipetman.

## D1200 Diamond Tips

- **New slender and long tip geometry** prevents the tip holder from coming in contact with wells and tubes.
- Designed for use even in long and narrow tubes.
- **Graduations** at 100  $\mu$ L and 1000  $\mu$ L make it easy to visually check the volume pipetted.



Specially designed for Gilson Pipetman Concept C1200, Pipetman Concept 8x1200 and 12x1200. Compatible with Pipetman P1000 and Pipetman Ultra U1000.